

A Safer, Smarter Way to Perform Infant Car Seat Testing

NOW AVAILABLE: Hospital discharge equipment designed for accurate, ergonomic infant car seat testing.



Nurses frequently end up kneeling on the floor with rolled up towels or improvising with beds and countertops.

The Infant Car Seat Tester from Patients First Products provides a safe, standardized platform that improves accuracy, infection control, and staff ergonomics.



Why Nurses Love the Car Seat Tester

Improves Testing Accuracy

Holds the infant carrier at the proper incline

Enhances Safety & Infection Control

Keeps carriers off floors, beds, and countertops

Reduces Strain on Staff

Eliminates kneeling, bending, & awkward positioning

Standardizes Discharge Protocols

Creates a repeatable, reliable testing process

Patients First Products: [Safety from the Start.](#)
See specs, videos, and more.

[Learn More at PFPCares.com](#)

SCAN FOR DETAILS



How It Works

The tester securely docks the infant carrier and allows easy positioning at the correct angle. Clinical staff can monitor the infant comfortably and consistently throughout the test.

Key Features



● Adjustable Testing Angle

- Adjustable incline for proper positioning
- Magnetic angle gauge for quick verification



● Ergonomic for Staff

- Positions carrier at a comfortable working height
- Eliminates floor and countertop testing



● Mobile & Stable Design

- Portable base for easy movement between rooms
- Locking brakes for stability during testing

● Fits All Standard Infant Carriers

- Secure straps hold the carrier in place
- Stable platform for testing procedures

● Built for Healthcare

- Designed for hospital environments
- Easy-to-clean materials

Improve Patient Experience with Safety from the Start

REQUEST A QUOTE TODAY

PHONE
239.220.0347

EMAIL
Info@pfpcares.com